## **PierenPlus basic**

Gebrauchsanweisung  
Operating instructions  
Mode d'emploi  
Gebruiksaanwijzing  
Használati utasítás

## PierenPlus basic

<b>D</b>	Gebrauchsanweisung.....	4
<b>GB</b>	Operating instructions.....	42
<b>F</b>	Mode d'emploi .....	68
<b>NL</b>	Gebruiksaanwijzing.....	94
<b>H</b>	Használati utasítás .....	120

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## Contents

Product features.....	43
Intended use.....	43
To ensure safe usage of the product.....	43
Operating elements.....	45
Keyboard.....	46
Replacing batteries.....	46
Operating device.....	47
Therapy mode.....	48
Programmes.....	49
Treatment examples.....	53
Therapy time mode.....	60
Date/time mode.....	61
Readout mode.....	62
Care and cleaning.....	62
Classification.....	63
New settings, alterations and repairs.....	63
Guarantee.....	63
Combination.....	63
Scope of supply.....	64
Accessories.....	64
Symbols.....	31
Technical data.....	31

## Product features

1. 15 programmes
2. Large LCD display for clear illustration of all functions and parameters:
  - Current programme
  - Duration of therapy
  - Intensity (two separate channels)
  - Frequency
  - Pulse width
3. Automatic switch-off when the PierenPlus basic is not in use.
4. Simple setting of intensity via keyboard.
5. Battery level recognition (the battery symbol appears where battery charging level is less than 3.7 V).
6. Electrode recognition  
(the intensity is reset to zero where no load is recognised).
7. Locking of current programme
8. Survey and readout of different operating parameters.

## Intended use

The PierenPlus basic serves the transcutaneous electrical stimulation of nerves and muscles in humans. Do not use the PierenPlus basic for any other purpose.

## To ensure safe usage of the product

Please read the instructions for use carefully prior to using the product!

- The PierenPlus basic may only be used with original accessories.
- Keep PierenPlus basic away from water and other liquids.

- Do not allow PierenPlus basic to be dropped, handled incorrectly, exposed to extreme temperatures or high humidity (only use at temperatures between 10 °C and 40 °C and a relative humidity of less than 90 %).
- Never use the PierenPlus basic if it functions incorrectly or has been damaged.
- Care must be taken when the PierenPlus basic is used on or in the proximity of children.
- Store the PierenPlus basic in the original packaging after use to protect it against damage and soiling.

### Precautionary measures

In the case of patients with an implanted electronic device (e. g. pacemaker) no stimulation should be carried out with the PierenPlus basic without the prior consultation of a doctor. The PierenPlus basic may only be connected to one patient for stimulation.

### Warning!

The simultaneous connection of the patient to a high-frequency surgical device may lead to burns beneath the electrodes. Operation in the vicinity of a short wave or microwave device may cause fluctuations in the output value of the PierenPlus basic. Portable and mobile RF communication devices may also influence the PierenPlus basic. The minimum surface of the electrodes should not be less than 2 cm<sup>2</sup>.

## Operating elements

The PierenPlus basic was designed for the stimulation of nerves and muscles in humans. All settings can be performed via the keys. The various operational status levels are illustrated in the display.

## Keyboard

The PierenPlus basic has seven keys:

- [ • ] Key for switching PierenPlus basic on and off
- [ P ] Key for programme selection
- [ E ] Key for recalling operating parameters
- [ ▲ ] Key for increasing intensity (left for Channel 1, right for Channel 2)
- [ ▼ ] Key for reducing intensity (left for Channel 1, right for Channel 2)

**Note: It is not possible to increase the intensity of both channels simultaneously.**

## Replacing batteries

The voltage of the batteries is monitored by the PierenPlus basic. If the voltage sinks below 3.7 V, then the battery symbol [ 🔋 ] appears in the display. It is now necessary to place four new 1.5 V batteries (AAA e.g. LR03) in the device.

- Switch off the PierenPlus basic.
- Open the cover of the battery compartment on the rear by sliding in the direction of the arrow.
- Remove the four used batteries.
- Insert four new batteries into the battery compartment.  
Please use the diagram in the battery compartment as orientation.
- Close the battery compartment again with the battery compartment cover.

Always remove the batteries if the device is not to be used for long periods. If rechargeable batteries are used, the instructions of the battery charger are to be followed.

Due to battery ordinance guidelines you are obliged to return used batteries or rechargeable batteries to your dealer or municipal collection point.

## Operating the device

### Start up

Prior to starting up the product the batteries supplied should be placed in the battery compartment (see also "Replacing batteries" section). All symbols appear briefly on the LCD display and a signal tone is heard.

### Switch on

Push the [●] key: You hear an acoustic signal and the PierenPlus basic is in standby mode.

### Ready mode

The LCD display shows the selected programme, the duration of therapy, the frequency, the pulse width and intensity.

1. Press the [ P ] key for programme selection.
2. Press the left [▲] key: the PierenPlus basic goes into therapy mode and the intensity for Channel 1 is increased by 1 mA (see Therapy-Mode).
3. Press the right [▲] key: the PierenPlus basic goes into therapy mode and the intensity for Channel 2 is increased by 1 mA.

### ATTENTION!

The device has a so called „electrode control". At an intensity of 10 mA, the intensity is put down to 0 if the electrodes are not connected properly.

4. Pushing and holding the left [▼] and the [ P ] button for three seconds:  
The [ P ] button is locked respectively unlocked. If the device is locked, the [ ] symbol appears in the display, an acoustic signal will sound up and no other programme can be selected.
5. Pushing the [ E ] button sets the device into therapy time mode (see therapy time mode). Pushing and holding the [ E ] button for three seconds leads to the date/time mode (see date/time mode).
6. Pushing and holding the left [▼] and [ E ] button for three seconds sets the device into readout mode (see readout mode).

**NOTE:** The removal of the battery resets the stimulation time to zero.

### Turn-off

Press [●] to turn-off the device: A single beep indicates that the PierenPlus basic turns-off. The PierenPlus basic turns-off automatically after 2 minutes if no key is pressed.

## Therapy mode

In therapy mode, the LCD display shows the selected programme, the remaining treatment time, the frequency, the pulse width and the intensity.

1. Push [●] or [ P ] to stop the therapy. The PierenPlus basic returns to standby-by mode
2. Push the left [▲] respectively the left [▼] button to adjust the intensity of channel 1 (from 0 to 60 mA).
3. Push the right [▲] respectively the right [▼] button to adjust the intensity of channel 2 (from 0 to 60 mA).
4. When the therapy time is over the therapy is interrupted and the device returns to ready mode.

## Programmes

The 12 fixed programmes possess the following parameters: Within this programme the therapy time can be adjusted from 10 to 90 minutes.

### Programme 1: Gate control 1

Both channels are operated simultaneously with the same frequency and pulse width.

- Frequency: 100 Hz
- Pulse width: 200  $\mu$ s
- Duration of therapy: 30 min

### Programme 1: Gate control 2

Both channels are operated simultaneously with the same frequency and pulse width.

- Frequency: 80 Hz
- Pulse width: 150  $\mu$ s
- Duration of therapy: 30 min

### Programme 3: Low frequency

Both channels are operated simultaneously with the same frequency and pulse width.

- Frequency: 2 Hz
- Pulse width: 250  $\mu$ s
- Duration of therapy: 30 min

#### Programme 4: Gate control + low frequency

The two channels work with different frequencies but the same pulse width (duration of therapy: 30 min).

##### Channel 1:

- Frequency: 100 Hz
- Pulse width: 200  $\mu$ s

##### Channel 2:

- Frequency: 2 Hz
- Pulse width: 200  $\mu$ s

#### Programme 5: Adaptation

Following completion of phase 1 there is an automatic transition to phase 2 (duration of therapy: 30 min).

##### Phase 1:

- Frequency: 100 Hz
- Pulse width: 150  $\mu$ s
- Working time: 10 min

##### Phase 2:

- Frequency: 2 Hz
- Pulse width: 200  $\mu$ s
- Working time: 20 min

#### Programme 6: Han

The device alternates continuously between phase 1 and phase 2 (duration of therapy: 30 min).

##### Phase 1:

- Frequency: 100 Hz
- Pulse width: 150  $\mu$ s
- Working time: 3 s

##### Phase 2:

- Frequency: 2 Hz
- Pulse width: 200  $\mu$ s
- Working time: 3 s

#### Programme 7: Burst

Pulse packages (bursts) are generated with the specified parameters.

- Frequency: 100 Hz
- Pulse width: 150  $\mu$ s
- Modulation: Burst 0.5 s (0.25 s + 0.25 s)
- Duration of therapy: 30 min

### Programme 8: Modulation

The frequency is changed continuously within the stated parameter borders.

The pulse width is adapted accordingly.

- Minimal frequency: 2 Hz  
→ Pulse width: 200  $\mu$ s
- Maximal frequency: 80 Hz  
→ Pulse width: 100  $\mu$ s
- Duration of modulation: 7.5 s
- Duration of therapy: 30 min

### Programme 9: Muscle training

The intensity is altered continuously in accordance with the stated parameters.

- Frequency: 50 Hz
- Pulse width: 250  $\mu$ s
- Intensity modulation:  
2 s rise → 5 s working → 1 s fall → 12 s pause
- Duration of therapy: 30 min

### Programme 10: Gate control dynamic 1

The intensity is altered continuously in accordance with the stated parameter, whereby both channels operate at delayed times to one another (massage).

- Frequency: 80 Hz
- Pulse width: 150  $\mu$ s
- Intensity modulation: 1 s rise → 1 s fall
- Duration of therapy: 30 min

### Programme 11: Gate control dynamic 2

The intensity is altered continuously in accordance with the stated parameter, whereby both channels operate at delayed times to one another (massage).

- Frequency: 80 Hz
- Pulse width: 150  $\mu$ s
- Intensity modulation: 0.25 s rise  $\rightarrow$  0.25 s fall
- Duration of therapy: 30 min

### Programme 12: Deep TENS

The programme works as in Programme 1. However, instead of one pulse there are four short pulses generated with the stated parameters.

- Frequency: 100 Hz
- Pulse width: 75  $\mu$ s
- Modulation: four pulses in a group, pulse width 75  $\mu$ s, duration of periods 275  $\mu$ s
- Duration of therapy: 30 min

## Programmes 13-15 are individually editable ("User-Programmes")

### Programme 13: User-Programme "Gate Control 1"

Both channels have continuously the same frequency and pulse width.

- Frequency from 1-120 Hz
- Pulse width from 75-300  $\mu$ s
- Therapy time from 10-90 minutes

### Programme 14: User-Programme "Muscle Training"

The intensity continuously changes according to the indicated parameters.

- Frequency from 1-120 Hz
- Pulse width from 75-300  $\mu$ s
- Intensity modulation:  
2 s rise  $\rightarrow$  5 s working phase  $\rightarrow$  1 s fall  $\rightarrow$  12 s pause
- Therapy time from 10-90 minutes

### Programme 15:

User-Programme "Gate Control dynamic 1" The intensity continuously changes according to the indicated parameters. The two channels work delayed (massage effect).


- Frequency from 1-120 Hz
- Pulse width from 75-300  $\mu$ s
- Intensity modulation: 1 s rise  $\rightarrow$  1 s fall
- Therapy time from 10-90 minutes

## Treatment examples

All electrode applications displayed and described are examples only and intended for initial orientation. The precise electrode application must be tested on a case-by-case basis.

### Tension headaches

Fig. 1: Programme 1, 2, 7

### Programme 1, 2, 5, 7

Tension headaches often occur on both sides and are often perceived as a painful hood, extending over the head from the nape of the neck. Here it is advisable to apply the electrodes to the cervical vertebrae muscle as shown

Fig. 2: Programme 5, 7


Fig. 3: Programme 2


Fig. 4: Programme 5, 7



Fig. 5: Programme 5, 7

in fig.1, preferably at points that react painfully to pressure. As a preventative measure, regular treatment of the shoulder-neck area may be helpful. The electrode application for this is illustrated in fig. 2.

## Migraine

### Programme 2


The migraine is a headache form usually occurring on one side, often accompanied by nausea, sensitivity to light and other symptoms. Clear points of pain often occur, such as in the region of the eyes or on the temples. During stimulation one electrode is applied directly onto the point of pain, above the eyebrow for pain in the vicinity of the eye. The second electrode can be placed on a second pain point or attached to the cervical vertebrae muscle on the same side (see fig. 3).

## Pain in the cervical vertebrae

### Programme 5, 7

Pain in the cervical vertebrae is often accompanied by limited mobility of the head and tense muscles. The pain tends to radiate from the shoulder area. During treatment with the PierenPlus basic the electrodes are attached to the cervical vertebrae musculature (see fig. 4). If the entire cervical vertebrae are inflamed or if pain radiates from both shoulders, then the two channels of the PierenPlus basic can be used to treat simultaneously, as illustrated in fig. 4 and 5.

Fig. 6: Programme 1, 3, 7


## Shoulder pain

Programme 1, 3, 7

Shoulder pain may emanate from the shoulder joint itself or from the surrounding structure. Treatment involves placing one electrode above the front armpit crease, the other above the rear armpit crease (see fig. 6).

Fig. 7: Programme 1, 7


## Shoulder/arm pain

Programme 1, 7

If the pain radiates from the shoulder down into the arm, then one electrode is applied to the main pain point on the shoulder and the other to the most painful point on the arm or the point on the arm where the pain ceases (see fig. 7).

Fig. 8: Programme 1, 7

## Elbow pain


Programme 1, 7

This is often the result of overloading syndrome. The outer or inner elbow is affected. During the treatment the electrodes are applied to the inner and outer ends of the elbow crease (see fig. 8) or, where the pain radiates into the forearm, one electrode on the main pain point of the elbow and a second on the radiating point of the forearm (see fig. 9).

Fig. 9: Programme 1, 7

Fig. 10: Programme 5


## Pain in the lumbar region, small of the back, sciatica

### Programme 1, 2, 5


Pain in the lumbar region is often the result of wear on vertebrae and spinal discs. This may lead to pressure on protruding nerves and radiating pain, e.g. in a leg. In the treatment the electrodes are attached to both sides of the vertebrae, onto the painful and tense muscle (see fig. 10). If the area of pain is extended, then both channels of the PierenPlus basic may be used (see fig. 11). If pain is radiated into the leg, then one electrode is applied to the main pain point in the back and the other to the main point of radiated pain or to the point on the leg where the radiated pain ends (see fig. 12).

Fig. 11: Programme 1, 2, 5

This electrode placement is particularly effective for sciatic pain.

Fig. 12: Programme 1, 2

## Hip pain

### Programme 1, 2, 3, 7

Pain in the hip area is often the result of wearing of the joints (arthrosis). Pain in the hip region may be treated with the application of electrodes above and below the joint (see fig. 13). In addition, pain is also radiated outwards, e.g. towards the groin. In this case an electrode is applied to the main point of pain on the hip and a second in the direction that the pain is radiated (see fig. 14).

Fig. 13: Programme 1, 3


Fig. 14: Programme 2, 7


Fig. 15: Programme 1, 2


Fig. 16: Programme 1



Fig. 17: Programme 1

## Knee joint problems

### Programme 1, 2

Common causes of problems in the knee joints are damage to the meniscus or the articular cartilage. One proven electrode application is to place the electrodes on the interior and exterior sides of the knee (see fig. 15). If there is a particular point of pain, such as the rear knee joint crease, then one of the electrodes may also be applied there. If both knees are painful, they may be treated simultaneously using the two channels of the PierenPlus basic.

## Trigeminal neuralgia

### Programme 1

Trigeminal neuralgia occurs in bursts of extremely unpleasant, shooting pain. Usually only one of the three trigeminal branches is affected. One electrode is applied directly in front of the ear, the second (depending on which branch is affected) above the eye (1st branch), beneath the eye at cheekbone height (2nd branch) or, also on the line of the eye, on the lower jaw (3rd branch) (see fig. 16). If the application of the electrodes or the stimulation itself is highly unpleasant, the healthy side of the patient is first treated with the same electrode application position as that foreseen for the side with the disorder. The TENS-specific effects also generate a pain-relieving effect on the side affected, even though this is not being directly treated!

Fig. 18: Programme 1


Fig. 19: Programme 9


Fig. 20: Programme 9



Fig. 21: Programme 9

## Pain following a shingles infection (post-zoster neuralgia)

### Programme 1

Significant pain may remain in the upper body and facial area following recovery from shingles. A high-frequency TENS treatment may have a significant effect on reducing this pain. One electrode is placed directly alongside the spinal column in the area affected, with the other electrode positioned at the limit of the pain radiation in the segment concerned, often on the stomach side (see fig. 17). If the application of the electrodes or the stimulation itself is highly unpleasant, the healthy side of the patient is first treated with the same electrode application position as that foreseen for the side with the disorder. The TENS-specific effects also generate a pain-relieving effect on the side affected, even though this is not being directly treated!

## Encouraging blood circulation

### Programme 3

This stimulation current application on the hand (see fig. 18) has a positive effect on blood circulation around the whole body! A stimulation of at least 20 minutes is required. It is possible to stimulate for up to 60 minutes without risk.

## Muscle stimulation in general

### Programme 9

Programme 9 contains four different modulation

Fig. 22: Programme 9


Fig. 23: Programme 9


Fig. 24: Programme 9



Fig. 25: Programme 9

phases (see programme overview). Modulation is the course of a stimulation pulse over a period of time.

A distinction is made between:

- rise time  
(duration until maximum current flow is achieved)
- working time (maximum current flow is available)
- fall time (duration until no more current flows)
- and pause time (no current flow).


With this breakdown into different times it is possible to adapt optimally to different physical conditions, that is, the natural processes in the muscle are stimulated. The rise time serves as a preparation for the musculature in the sense of a pre-tensing. During the working time the tension of the muscle is maintained at a constant level. The pause time is particularly important. It is common knowledge that muscles need time to recover from strenuous use. This must also be taken into consideration during electro stimulation. Following application of the electrodes the device is switched on and the current strength slowly regulated upwards, with a tangible muscular response, without this being uncomfortable or painful.

Muscle stimulation under medical conditions:

- Inactivity atrophy following long periods of immobility of limbs or other parts of the body
- Muscular atrophy due to myopathy

Abb. 26: Programme 9

- Muscular imbalance and insufficiency due to orthopaedic reasons, e.g. myostatic insufficiency, chronic chondropathia patellae
  - Postoperative functional dysfunction, e.g. following the fitting of an artificial hip
  - Trauma-related functional dysfunction
  - Functional dysfunction following cerebral paralysis
  - Stimulation of facial musculature with facial paresis
- Preoperative muscle training in anticipation of dysfunction
  - Intention training in the case of co-ordination dysfunction following peripheral palsy or other neurological disorders
  - Build-up of innervated muscle fibre following peripheral paralysis
  - Urinary or stool incontinence in the case of pelvic floor insufficiency and innervated musculature
  - Scoliosis for the targeted stimulation of convex back musculature These are highly specific applications which should only be performed following consultation with a doctor or physiotherapist. The correct placement of the electrodes is of particular importance here. The stimulation duration must also be individually specified.

### General rules are:

When beginning muscle stimulation therapy do not work longer than 15 minutes, as this may lead to muscle ache and a deterioration in the complaint. Over a period of 3 to 4 days the period can be gradually increased to up to 60 minutes.

## Therapy time mode

In therapy time mode, the LCD display shows the selected programme, the frequency, the pulse width and the intensity. The therapy time flashes. Pushing the

left or right [▲] button increases the therapy time by 5 minutes (max. 90 minutes). Pushing the left or right [▼] button decreases the therapy time by 5 minutes (min. 10 minutes). If the [ E ] or [●] button is pushed in the programme 1-12, the set time is stored and the therapy time mode is left.

Pushing the [ E ] or [●] button in programme 13-15 puts the device in set-up mode for the therapy frequency. The therapy frequency flashes. Pushing the left [▲] button increases the frequency by one Hertz (max. 120 Hz). Pushing the left [▼] button decreases the frequency by one Hertz (min. 1 Hz).

Pushing the [ E ] button in programme 13-15 puts the device in set-up mode for the pulse width. The pulse width flashed. Pushing the left [▲] button increases the pulse width by five microseconds (max. 300  $\mu$ s). Pushing the left [▼] button decreases the pulse width by five microseconds (min. 75  $\mu$ s). Pushing the [ E ] or [●] button will store the set parameters and the therapy time mode is left.

The PierenPlus basic turns-off automatically after two minutes if no key is pressed. In this case no parameters will be stored.

## Date/time mode


The LCD display shows two figures. The left figure indicates the date, the right figure indicates the time. Pushing the left [▲] button increases the date by one day (max. 30). Pushing the left [▼] button decreases the date by one day (min. 1). Pushing the right [▲] button increases the time by one hour (max. 23). Pushing the right [▼] button decreases the time by one hour (min. 0). By pushing the [●] button all parameters are stored and the date/time mode is left.

The PierenPlus basic turns-off automatically after two minutes if no key is pressed. In this case no parameters will be stored.

## Readout mode

If no values have been stored at the data memory of the PierenPlus basic, the LCD display will show „NULL“ for two seconds and return to standby mode. Otherwise the LCD display shows the parameters used during the last therapy session (see fig.).

- Selected programme = P-01
- Frequency = 100 Hz
- Pulse width = 200 µs
- Therapy time = 0 min 44 s
- Therapy day = 5.
- Therapy session = 3.

You have the possibility to store 90 parameters, 30 days and 3 parameters each day. Pushing the left [▼] button shows the operating parameters of the last therapy day (until the first). Pushing the left [▲] button shows the operating parameters of the next day (up to the latest). Pushing the right [▼] button shows the operating parameters of the last therapy session (until the first of the respective day). Pushing the right [▲] button shows the operating parameters of the next therapy session (until the last of the respective day). Pushing and holding the [E] for five minutes deletes the data memory. The LCD display shows „CLr“ for two seconds and the PierenPlus basic returns to standby mode.

The PierenPlus basic turns-off automatically after two minutes if no key is pressed.

## Care and cleaning

No special care and cleaning materials are required for the PierenPlus basic. The PierenPlus basic can be cleaned with a soft, lint-free cloth. It is vital to ensure

that no moisture enters into the interior of the device. Should moisture penetrate despite this, the device must be subjected to a technical inspection prior to being used.

## **Classification**

The PierenPlus basic is assessed as Class IIa in accordance with Appendix IX of the EU directive on medical products.

## **New settings, alterations and repairs**

The manufacturer only holds itself responsible for the safety and performance of the PierenPlus basic when new settings, alterations and repairs are carried out by a person authorised by the manufacturer and the PierenPlus basic is operated in accordance with the operating instructions.

## **Guarantee**

We issue a guarantee of one year on the PierenPlus basic, with the exception of expendable items such as cable, electrodes etc.

## **Combination**

The PierenPlus basic may be combined with all articles included in the scope of supply and accessories.

## Scope of supply

Quantity	Article	REF
1	PierenPlus basic	462061
1	Cable Type 7 (VPE = 2 units)	104741
1	Self-adhesive electrode 50 mm x 50 mm (VPE = 4 units)	281006
4	1.5 V battery LR03 micro (AAA)	450780-0011
1	Operating instructions	451600-0241

## Accessories

Quantity	Article	REF
1	Self-adhesive electrode 32 mm round (VPE = 4 units)	281004
1	Self-adhesive electrode 50 mm round (VPE = 4 units)	281005
1	Self-adhesive electrode 50 mm x 50 mm (VPE = 4 units)	281006
1	Self-adhesive electrode 50 mm x 90 mm (VPE = 4 units)	281007
1	Self-adhesive electrode 50 mm x 130 mm (VPE = 4 units)	281008

## Symbols

Attention: Read Operating instructions!

BF type application part. Protection against electric shock.

The year of construction of the product follows this symbol.

**REF**

The article or order number of the product follows this symbol.

**SN**

This symbol is followed by the serial number of the device as well as the production run indicated by a letter.

This equipment is marked with the recycling symbol. It means that at the end of the life of the equipment you must dispose of it separately at an appropriate collection point and not place it in the normal domestic unsorted waste stream. This will benefit the environment for all.

**CE 0197**

Conform with Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

## Technical data

2-channel electro stimulator with separated outputs,  
12 integrated and 3 editable programmes.

Power supply:	6.0 V (4 x 1.5 V battery AAA type micro e. g. LR03) or 4.8 V (4 x 1.2 V rechargeable battery AAA)
Current consumption:	Max. 100 mA (therapy) Approx. 30 $\mu$ A (power down)
Dimensions:	Approx. 140 mm x 64 mm x 28 mm
Weight:	Approx. 96 g (without batteries)
Output current:	0-60 mA (with 1 k $\Omega$ load)
Pulse form:	Positive rectangle with negative component
Frequency range:	1-120 Hz
Pulse width:	75-300 $\mu$ s
Operating conditions:	Temperature range: 10 °C to 40 °C Relative humidity: 30 % to 90 %
Storage conditions:	Temperature range: -10 °C to 55 °C Relative humidity: 10 % to 90 %

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# PIERENKEMPER GmbH
Pierenkemper GmbH